

A Symposium on Cardiovascular Imaging

“New Horizon of Cardiac CT & MRI”

Seoul, October 13, 2006

SIEMENS

State-of-the Art Technology in Cardiac MR

Gerhard Laub, Ph.D.
Siemens Medical Solutions

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 1

Cardiovascular Tools

anatomy & morphology

coronary MRA

function & wall motion

perfusion

viability

angiography

Total Imaging Matrix: Tim

SIEMENS

- Up to 105 coil elements
- Connected into 32 channels

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 3

Parallel in all Directions

SIEMENS

- High-speed / high-resolution parallel imaging.
- In all directions, whole body.
- No need for specific PAT coils anymore.

Parallel imaging independent of slice orientation

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 4

Real-time Cine Imaging

SIEMENS

- TSENSE x 3
- 65 msec true temporal resolution
- 80x192 matrix, 255mm x 340mm x 8mm
- TE 1.1 msec, 80° flip angle TrueFISP

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 5

High spatial resolution in real time

32 channel cardiac coil, PAT 4, T-SENSE

SIEMENS

- Res: 1.1 x 0.8 mm²
- TA: 92 ms / image
- PAT 4

Courtesy of Rapid Biomedical, Wuerzburg Germany Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 6

A Symposium on Cardiovascular Imaging

“New Horizon of Cardiac CT & MRI”

Seoul, October 13, 2006

GRE-EPI with T-SENSE X 2 **SIEMENS**
 Case example: anterior and inferior first-pass defects

STRESS

REST

Courtesy of Northwestern Memorial Hospital, Chicago Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 7

3D Perfusion with iPAT2 and T-SENSE **SIEMENS**
 Full cardiac coverage with PAT 8 (4 x 2)

- GRE-EPI Inversion Recovery / Fat Saturation
- Matrix: 128x64; FA: 30°; BW/pixel: 1860 Hz; EPI 4; 1D T-SENSE R=2;
- 8 Partitions, 50% slice resolution; TI 200 ms; TR / TE 5.5/1.2ms

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 8

Two in one : Cine – Late Enhancement **SIEMENS**
 Combined wall motion and late enhancement in one exam

J. Kim, R. Sotser, R. White, A. Stillman, Cleveland Clinic Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 9

Cine Late Enhancement **SIEMENS**

Images courtesy Dr. Regenfus, Uni Erlangen Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 10

3D Coronary Imaging **SIEMENS**
 without contrast agent

- New navigator techniques
 - Adaptive motion correction
- New visualization tools

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 11

Free Breathing Coronary MRA **SIEMENS**

LAD Proximal RCA Distal RCA

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 12


A Symposium on Cardiovascular Imaging

“New Horizon of Cardiac CT & MRI”

Seoul, October 13, 2006

Free Breathing Coronary MRA

SIEMENS




Soap Bubble Reformat

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 13

Non-contrast MRA: Dilated aorta

SIEMENS

TR / TE / flip = 2.3 ms / 1.0 ms / 90°
 FOV = 500 mm × 500 mm
 Matrix = 320 × 320
 slic thk = 4 mm interpolated to 2 mm
 TA = 6min, 15sec
 Free breathing




Thickened aortic valve

Courtesy of Dr. Finn, UCLA, Los Angeles
 Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 14

Transposition of the great vessels

SIEMENS




Functional aorta Left ventricle
 Functional pulmonary artery Left ventricle
 Right ventricle

TR / TE / flip angle = 2.3 ms / 1.0 ms / 90°
 FOV = 420 mm × 420 mm, Matrix = 256 × 256
 Slice thickness = 4 mm interpolated to 2 mm
 TA = 9min, 59s.

Courtesy of Dr. Finn, UCLA, Los Angeles
 Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 15

Thorax-Abdomen

SIEMENS




Thin MIP MIP

Tim:
 Large FOV
 &
 High Resolution

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 16

Dynamic MRA with T-SENSE 2

SIEMENS



Acceleration: **8**
 Matrix 256x256x40
 FOV 400x400x120 mm³
 TR/TE 3.5/1.1 ms
 TSENSE 4x2
 Acq time/meas **2.7 s**
 Rapid 32 channel coil
 Tim Avanto [76x32]


Courtesy of P Kellman, NHLBI, NIH, Bethesda, MD, USA
 Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 17

Dynamic MRA Using TWIST:

SIEMENS

As is well-known, k-space can be divided into two regions, A and B. “A” defines the overall image contrast, and “B” adds object details.

In the TWIST sequence, region B is scanned at a lower sampling rate to increase the frame rate.




Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 18

A Symposium on Cardiovascular Imaging

“New Horizon of Cardiac CT & MRI”

Seoul, October 13, 2006

Dynamic MRA Using TWIST: SIEMENS



Congenital Vascular Disease

Time-resolved (iPAT x 2)
1.3 x 1.1 x 4.0mm³
1.8 sec / frame

Avanto, UCLA Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 19


Dynamic Peripheral MRA at 3.0 T SIEMENS



AVM, time-resolved MRA using TWIST

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 20

Cardiac MR at 3T SIEMENS



1.5T


3.0T

1.0 x 1.0 x 3 mm

Courtesy of Dr. Miller, University of Tuebingen Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 21

Comparison of 3T and 1.5T SIEMENS

- resting first-pass perfusion -



1.5T


3.0T

→ Higher SNR at 3T

Courtesy of Drs. J. Salanitri and J. Carr, Northwestern University Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 22

Comparison of 3T and 1.5T SIEMENS

- tagging sequence -



1.5 T


3.0 T

TT = 145 ms TT = 425 ms TT = 700 ms

→ Better CNR at 3T

Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 23

Cardiac Function @ 3T SIEMENS



Cine TrueFISP


courtesy of Dr. Miller, University of Tuebingen Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 24

A Symposium on Cardiovascular Imaging

"New Horizon of Cardiac CT & MRI"

Seoul, October 13, 2006

Cardiac MRI @ 3T SIEMENS




Aortic valve leaflet defect Aortic regurgitation

Cine FLASH to avoid off-resonance artifacts

Courtesy of Dr. Finn, UCLA Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 25

Viability & Function @ 3T SIEMENS




Cine TrueFISP IR single shot TurboFLASH

IR segmented TurboFLASH (4 hb) IR segmented TurboFLASH (4 hb)

courtesy of Dr. DiBello, University of Utah Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 26

Cine, Perfusion, and Viability @ 3T SIEMENS




Cine TrueFISP SR TurboFLASH + TSENSE Phase sensitive IR TurboFLASH

function perfusion viability

Courtesy of Dr. Finn, UCLA Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 27


Pulmonary MRA SIEMENS



high res, 20sec Dynamic MRA (TWIST) 1sec resolution

Courtesy of Dr. Finn, UCLA Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 28


Aortic Dissection SIEMENS



Dynamic MRA (TWIST) High-res MRA (Thin MIP)

Courtesy of Dr. Finn, UCLA Copyright © 2006 Siemens Medical Solutions USA, Inc. All rights reserved. 29

State-of-the Art Cardiovascular Tools 1.5T and 3T



anatomy & morphology dynamic MRA

function & wall motion perfusion viability angiography