


Coronary MR Angiography

Hajime Sakuma, MD

Department of Diagnostic Radiology
Mie University Hospital
2-174 Edobashi, Tsu, Mie 514-8507, Japan
sakuma @clin.medic.mie-u.ac.jp

Advantages of coronary MRA

- No ionizing radiation
- No contrast injection
- No interference by calcification



Weakness of MRI


- Substantially slower than 64-slice MDCT
 - Acquisition time of 3D volume (FOV 32x32x20cm, acquisition resolution 0.6x0.6x0.6mm, no ECG gating)
 - 64-slice MDCT < 2 sec.
 - MRI > 2 min.
- 40,960 phase encoding steps
TR 3ms, SENSE factor 2, half scan

MR detection of coronary artery stenoses


Investigator	Reference	Method	Sensitivity (%)	Specificity (%)
Manning WJ	NEJM 1993;328:828	2D breath-hold	90	92
Duerinckx AJ	Radiology 1994;193:731	↑	63	82
Van Geuns RJ	Radiology 2000;217:270	3D breath-hold	68	97
Regenfus M	JACC 2000;36:44	↑	94	57
Woodard PK	AJR 1998;170:883	3D navigator, retrospective	73	N/A
Sandstede JJ	AJR 1999;172:135	↑	81	89
Sardanelli F	Radiology 2000;214:808	↑	82	89
Kim WY	NEJM 2001;345:1863	3D navigator, prospective	93	42
Bogaert J	Radiology 2003;226:770	↑	50	90
Kefer J	JACC 2005;46:92	↑	75	77

1993-2005 Sensitivity 50-94%, Specificity 42-97%


Breath-hold 2D coronary MRA



Edelman RR, Radiology 1991;181: 641-643
Manning WJ, N Eng J Med 1993;328: 828-832
Sakuma H. AJR 1994; 163:533-537




Respiratory gated 3D coronary MRA



Kim WJ, New Engl J Med 2001; 345:1863
Sardanelli F. Radiology 2000;214:808
Sakuma H. Invest Radiol. 1999;34:503




Breath-hold 3D coronary MRA



3D Breath-hold Coronary MRA - Multi-center study -

Foo TK, Ho V, Saranathan M, Cheng L, Sakuma H, et al.


- Breath-hold 3D FIESTA
- 1.1x1.2x2 mm resolution
- 77%(79/103) coronary arteries visualized
- In 79 arteries visualized, 82% showed diagnostic image quality (Grade 3-4)




GE CVI, Mie Univ.

Radiology 2005;235:1025-1030

Whole Heart Coronary MRA



Whole heart coronary MR angiography using steady state free precession




5:36 min/volume x 4 (Targeted volume)

13:48 min x 1 (Whole heart)


Weber OM. *MRM* 2003;50:1223
Sakuma H, *Radiology* 2005;237:316

Free-breathing 3D coronary MRA

- Navigator echo
 - prospective real-time respiratory gating
 - gating window = +/- 2.5mm
 - adoptive motion correction



Whole heart coronary MRA in 131 patients using optimized acquisition window in the cardiac cycle



During diastole
83/131 (63.3%) patients
HR = 65.1 BPM ± 9.4
Window = 152 ms ± 67


During systole
48/131 (36.7%) patients
HR = 83.2 BPM ± 9.5
Window = 98 ms ± 26

Axial cine MRI (50 phases)


Sakuma H. *J Am Coll Cardiol* 2006 (in press)

Whole heart coronary MRA

Sequence	3D B-TFE
FOV	28x28x12cm
Preparation	T2 prep SPIR
Acq. Mat.	256x256x80
TR	4.6ms
TE	2.3ms
Flip angle	90 deg
SENSE fact.	2
Imaging time	6-20 min




Whole heart coronary MRA Normal coronary artery



Volume Rendering

Whole heart coronary MRA 61M



Curved MPR of RCA

X-ray coronary angiography

Whole heart coronary MRA 61M




Initial results

- **Imaging time**
13.8min ± 3.8 (range 6.8-23.8 min)
- **Imaging completed** 34/39 (87.2%)
- **Study success rate** 32/39 (82.1%)
High quality (score ≥3 in 4 point-scale) images
in all three major coronary arteries


Sakuma H, Ichikawa Y, *Radiology* 2005;237:316

78-year-old man with LAD stenosis (HR = 64/min)



A and B, Curved multiplanar reconstruction of whole heart coronary MRA.
C and D, Volume rendering of whole heart coronary MRA.
E, Catheter X-ray coronary angiography

**54-year-old man with LCX stenosis
(HR = 78/min.)**



A, Curved multiplanar reconstruction of whole heart coronary MRA.
B, Volume rendering of whole heart coronary MRA.
C, Catheter X-ray coronary angiography.

Detection of significant coronary artery stenoses with whole heart coronary MRA using optimal acquisition window in the cardiac cycle in 113 patients

	n	Sensitivity, %	Specificity, %	PPV, %	NPV, %	Accuracy, %
Per patient	113	82(69-91)	90(79-96)	88(74-95)	86(75-93)	87(79-92)
Per vessel	452	78(66-86)	96(93-97)	79(67-87)	95(93-97)	93(90-95)
RCA	113	85(65-95)	95(88-98)	85(65-95)	95(88-98)	93(86-97)
LM	113	NA	98(93-100)	NA	100(96-100)	98(93-100)
LAD	113	77(56-90)	95(88-99)	83(62-95)	93(85-97)	91(84-95)
LCX	113	70(47-86)	93(86-97)	73(50-88)	92(84-97)	89(81-94)
Per segment	1000	78(68-85)	96(95-97)	69(60-77)	98(96-98)	94(93-96)


Study success rate = 86% (113/131)

All arteries with diameter of ≥ 2 mm on coronary angiography were evaluated.

Sakuma H, Ichikawa Y, et al. J Am Coll Cardiol 2006 (in press)

How we use coronary MRA in the era of 64-slice MDCT.

64-slice MDCT angiography



76M RCA occlusion

Light Speed VCT64, Virtual Place
Matsusaka Central Hospital, Mie, Japan

Coronary MRA vs. CTA

- Kefer J, JACC 2005;46:92
- Coronary MRA - free breathing, 3D balanced TFE on targeted double oblique slices

	Sensitivity	Specificity
Coronary MRA	75%	77%
16-slice CTA	82%*	79%*

* p=N.S.

Diagnostic accuracy of CTA and MRA in detecting patients with significant CAD

Whole heart coronary MRA

	Sensitivity	Specificity
Jahnke C, Eur Heart J 2005 ;26:2313	78%	91%
Sakuma H, JACC 2006 (in press)	82%	90%


64-slice MDCT

	Sensitivity	Specificity
Mollet NR. Circulation 2005;112:2318	100%	92%
Pugliese S. Eur Radiol 2005;16:1	100%	90%
Leber AW. JACC 2005;46:155	90%	95%

How we use coronary MRA in the era of 64-slice MDCT.

1. Anomalous coronary artery
2. Coronary aneurysm in Kawasaki disease.
3. Patients with renal failure.
4. Patients with heavy coronary calcification.
5. As a part of one-stop-shop CMR study
6. Screening CAD in asymptomatic subjects


31-year-old man, VT after exercise



Anomalous origin of RCA

Whole heart coronary MRA

- 71-year-old man, chest pain on effort
- Coronary CT angiography was performed using 64-slice MDCT




Volume rendering

Angiographic view

Curved MPR

Coronary CT Angiography (GE VCT64)

- 71-year-old man, chest pain on effort
- Whole heart coronary MR angiography was subsequently performed.




Volume rendering

Curved MPR

X-ray coronary angiography

Whole heart coronary MR angiography

MRA as a part of one-stop-shop cardiac MR study



Cine MRI

Stress-rest

Delayed enhanced


Whole heart 3D

LV function

Perfusion

Viability

Morphology



Acknowledgements

Mie University Hospital
 Department of Radiology
 Nanaka Ishida, MD
 Kan Takeda, MD
 Shinichi Takase, RT
 Genji Asanuma, RT
 Department of Cardiology
 Tairo Kurita, MD
 Katsuya Onishi, MD
 Masaaki Ito, MD

Matsusaka Central Hospital
 Department of Radiology
 Yasutaka Ichikawa, MD
 Shuji Chino, MD
 Tadanori Hirano, MD
 Chikai Shoji, RT
 Tetsuya Ochiai, RT
 Department of Cardiology
 Katsutoshi Makino, MD

Philips Medical Systems